


PPE holder

PRODUCT FACT SHEET

8



06/2024



MADE IN SWEDEN
www.silentia.com

Product specifications


The Silentia PPE (Personal Protective Equipment) holder gives easy access to hand sanitizers, gloves, wet wipes, and other products. The PPE holder mounts onto Silentia screens by hand without the need for tools and holds PPE accessories with Silentia accessorie clamps.

Material

The profiles are made of aluminum. The top and bottom connectors are made of polyamide reinforced with fiberglass. The sliding hook for PPE holder 1 is made of ABS.

Dimensions & Sizes


PPE holder 1, compatible with:
Bed end screen, mobile folding screen



Height of your screen


One size fits all.

Art.no 06580



No. of holders: _____

PPE holder 2, compatible with:
Folding screen



Height of your screen

1,45m ☐

Art.no 06584

1,65m ☐

Art.no 06585

1,85m ☐

Art.no 06586




max. load 1kg



No. of holders: _____

PPE holder 3, compatible with:
Semifix screen, fixed screen



Height of your screen

1,45m ☐

Art.no 06581

1,65m ☐

Art.no 06582

1,85m ☐

Art.no 06583



max. load 6kg



No. of holders: _____

Product specifications

Care & maintenance

EasyClean™ All products meet strict hygiene requirements and all surfaces are accessible for cleaning. The PPE holder can easily be removed from the screen and cleaned. The clamps can be cleaned separately in a dishwasher.

Handling/Use

The PPE holder easily clicks onto the Silentia screen. Slide the PPE holder along the screen for optimal positioning.

The PPE products are attached to the holder with accessory clamps, making it easy to swap or move products. The accessory clamps can easily slide along the holder to the preferred position.

Sustainability

Silentia focuses on developing long-lasting products that minimize our environmental impact. Since 1990, Silentia has been committed to improving patient care around the world with a safe, hygienic, and resource-efficient privacy system. [LINK TO SUSTAINABILITY PAGE](#)

Recycling & disassembly

The PPE holder is recyclable and all parts can be completely separated and individually delivered to a suitable recycling facility.

Certificates / Listing



EN ISO 9001:2015



EN ISO 14001:2015



Listed as interior solution for the
"WORLD'S GREENEST HOSPITAL"

Tests

Silentia can provide fire tests, durability tests, etc. Please contact the Silentia sales team for specific information.

Warranty

Silentia provides a 2-year warranty. Spare parts are available. Please contact the Silentia sales team for further information about warranty terms and conditions.

Downloads

[LINK TO 2D/3D FILES](#)

[LINK TO SILENTIA MANUAL](#)


[LINK TO CARE & MAINTENANCE VIDEO](#)

Installation

PPE accessories


Please contact the Silentia sales team for further information.

PPE accessory clamp: Enables PPE products to be attached to the holder and slide to the desired position.



Art.no 06590
No. of clamps: _____



PPE adapter plate: Enables any PPE product to be attached to the accessory clamp. (Requires accessory clamp.)



Art.no 06591
No. of plates:
(Mounting bolts included) _____


Attached load diagram

Please contact the Silentia sales team for further information.




Mounting

PPE holder 1, compatible with:
Bed end screen, mobile folding screen
max. load 4kg



PPE holder 2, compatible with:
Folding screen
max. load 1kg



PPE holder 3, compatible with:
Semifix screen, fixed screen
max. load 6kg

